SPEIKOCOPAL

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Lack: 12,5 % Copalharz, Diethylether, Ethanol

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Stand der Information: November 2005

Composition:

Varnish: 12.5 % copal resin, diethyl ether, ethanol

Applications:

For the lining of flat cavities.
For protection varnish of silicate and siliko-phosphate cements.

Contraindications:

The application of SPEIKOCOPAL is contraindicated in the presence of plastic fillings, glass-ionomer cements, and polycarboxylate cements, as well as in the presence of an allergy to any of the ingredients.

Side effects/Interactions:

No systemic side effects of resins and solvents when used as cavity varnishes are known to date. In the presence of certain plastic fillings, cavity varnishes may interfere with hardening.

Dosage instructions, form and duration of application:

After cleaning and drying the cavity, the varnish is removed from the storage vessel and introduced into the cavity. The varnish is applied to the cavity in a thin layer and dried in a gentle stream of air. If a second layer of the cavity varnish has to be applied, application must not take place until the first layer has dried completely.

Before the filling material is applied, any cavity varnish must be removed from the margin of the cavity.

Special instructions:

Diethyl ether is a very volatile, highly flammable liquid and is highly explosive in contact with oxygen and dinitrogen monoxide (nitrous oxide, laughing gas).

Diethyl ether vapour is heavier than air and therefore always collects at the lowest point. Radios, contacts and electrostatic charges as well as electrocautery may trigger an explosion. Under the influence of light and oxygen, diethyl ether may degrade explosively to toxic peroxides. It must therefore be stored in brown, light-impermeable bottles. As possible keep in a cool place.

This medical product must not be used after the expiry date.

Keep out of reach of children.

Available in the following packs:

SPEIKOCOPAL varnish Bottle 20 ml Article No. 1049

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SPEIKO® - Dr. Speier GmbH Postfach 7204 D-48039 Münster